

## Anesthesia Gas Exposure Monitor



OR personnel may be exposed to waste anesthetic gases escaping from various components of the delivery system. PACU personnel encounter exposure to these gases from out-breathing recovering patients. The Kem Medical Waste Anesthetic Monitor # 8640 measures Sevoflurane, Desflurane, and Isoflurane concentrations, and #8530 measures for Nitrous Oxide. Accurate, fully validated monitors are also lightweight and convenient to use. Price includes analysis and pre-paid mailer.

**Source URL (retrieved on 01/27/2015 - 1:27pm):**

[http://www.surgicalproductsmag.com/product-releases/2005/04/anesthesia-gas-exposure-monitor?qt-recent\\_blogs\\_articles=0](http://www.surgicalproductsmag.com/product-releases/2005/04/anesthesia-gas-exposure-monitor?qt-recent_blogs_articles=0)