


Post-Rhinoplasty Splint

Published on Surgical Products (<http://www.surgicalproductsmag.com>)

Post-Rhinoplasty Splint



The Jiffy Splint comes packaged in a clear pack with a colored sheath for easy removal. Once 13 cc of water is injected into the packet, it is inserted into the slot at the top of the Jiffy Heater. The splint is ready in three minutes, but may stay in heater for hours without consequence and still be ready at the exact moment the doctor needs it. The new adhesive layer and hole pattern allow for much longer wear and protection for the post-rhinoplasty patient.

Source URL (retrieved on 01/25/2015 - 8:36pm):

http://www.surgicalproductsmag.com/product-releases/2005/08/post-rhinoplasty-splint?qt-recent_videos=0&qt-digital_editions=0