

## Active Patient Warming System



Geratherm active warming systems are used for the prevention and treatment of hypothermia during patient transport or in the Perioperative setting. Five carbon fiber warming panels distribute radiant heat in a full size blanket or 2/3 sized blanket. The edges are heat sealed providing protection from bodily fluids and facilitates cleaning and disinfecting. The control unit has a warming temperature range of 98.6 F to 105 F with both audible and visual alarms and automatic system shut off for all faults detected included in the design.

**Source URL (retrieved on 01/29/2015 - 8:37pm):**

<http://www.surgicalproductsmag.com/product-releases/2007/09/active-patient-warming-system>