

Instrument Sterilization Bag



The new Time Out Bag is validated for use inside of instrument trays and containers to hold separate items, in compliance with AAMI ST79. In addition, printed clearly on the bag is the message “Time Out” alerting staff to take a critical pause during surgery preparation. This paper bag is made of 100 percent medical grade paper, which is porous to sterilant. Further, the moisture absorbent paper helps to dissipate condensate, helping to fight wet pack.

???

Source URL (retrieved on 01/26/2015 - 11:04pm):

http://www.surgicalproductsmag.com/product-releases/2008/09/instrument-sterilization-bag?qt-recent_blogs_articles=0&qt-digital_editions=0