

# Disposable Patient Repositioning Sheets



The Solo RepoSheet<sup>®</sup> is designed help control infection and cross-contamination. This disposable patient repositioning sheet is usually applied above the fitted sheet and below the draw sheet, but can be used in place of the draw sheet. Because the Solo RepoSheet rests under the patient ready to be used, the sheet is composed of 100 percent light weave breathable cotton for the patient's comfort. The sheet isintended to be used in combination with a facility's existing overhead (ceiling) or mobile floor lifts. In addition to eliminating rotational and compression forces associated with moving patients up in bed, the Solo RepoSheet can eliminate cumulative trauma associated with turning patients, holding them in side-lying, and rotating a patient from supine to prone.

**Source URL (retrieved on 01/28/2015 - 5:26am):**

[http://www.surgicalproductsmag.com/product-releases/2008/12/disposable-patient-repositioning-sheets?qt-recent\\_videos=0](http://www.surgicalproductsmag.com/product-releases/2008/12/disposable-patient-repositioning-sheets?qt-recent_videos=0)