

## Low-Friction Transfer & Repositioning Sheets



The Stat HandySheet<sup>®</sup>; and the Stat HandyTube<sup>®</sup>; help facilitate the placement of slings under patients or the movement of patients in bed. Manufactured of extremely low friction fabrics, these new products are designed to reduce skin abrasions to the patient while also minimizing the caregiver's workload during a variety of patient transfer and repositioning maneuvers. Stat HandySheets are designed to be used in pairs and are recommended to facilitate placement of a sling while a patient is sitting in a wheelchair or lying on a bed. The Stat HandyTube is designed to assist caregivers in repositioning, turning patients or performing lateral transfer maneuvers.

**Source URL (retrieved on 03/04/2015 - 10:00pm):**

<http://www.surgicalproductsmag.com/product-releases/2008/12/low-friction-transfer-repositioning-sheets>