

Reusable Return Electrode



The Mega Soft from Megadyne is set on the operating room table, eliminating the small disposable sticky pad that can damage the patient's skin. According to the company, the Mega Soft pad is large—extending at least half the length and full width of an adult patient's torso. The patient lies on the Mega Soft pad to achieve full and safe electrical patient contact. If the patient is in contact with only a small portion of the pad, safety still prevails because the Mega Soft is a self-contained current limiting device. Additional features include:

- A combination of the Mega Soft patient return electrode technology with an O.R. table pressure reduction pad made with Akton® Polymer, which are among the most effective in reducing direct pressure, shear, and friction which contribute to the formation of decubitus ulcers.
- The ability to be used in a variety of surgical procedures and positions providing a safe, easy to use return electrode combined with an effective O.R. table pressure reduction pad.

For more information, visit www.megadyne.com [1]

Source URL (retrieved on 02/01/2015 - 9:49am):

<http://www.surgicalproductsmag.com/product-releases/2010/04/reusable-return-electrode>

Links:

[1] <http://www.megadyne.com/>