

Flash-Guard



Millennium Surgical offers the Flash-Guard®, the only container system capable of sterilizing full instrument sets, 16 pounds or less, of non-porous, non-cannulated instruments at 270°F for 3 minutes. For cannulated and porous items, Flash-Guard meets the 10-minute cycle requirement, as well. These sterilization parameters are the recommended practices of AAMI and AORN.

For more information, visit www.surgicalinstruments.com [1].

Source URL (retrieved on 01/29/2015 - 1:44am):

http://www.surgicalproductsmag.com/product-releases/2012/11/flash-guard?qt-most_popular=0&qt-recent_blogs_articles=0&qt-recent_content=0

Links:

[1] <http://www.surgicalinstruments.com/home.asp?uri=1000&frst=1000>