

Laparoscopic Instrument Cannula Cleaning Device



InstruSponge[®] is designed to clean internal channels of take-apart laparoscopic instruments, and is attached to a flexible plastic wand that allows for easy maneuvering through complex internal channels, cannulas, and lumens to loosen and expel contaminants. InstruSponge's foam tip is impregnated with Endozime[®] AW Triple Plus[®] with APA (Advanced Proteolytic Action) to gently wipe clean the inside surface of the take-apart channel while instantly and actively digesting bioburden upon contact. InstruSponge's unique foam tip will not scar instrument walls, unlike conventional channel wire bristle brushes.

Source URL (retrieved on 01/28/2015 - 4:29pm):

http://www.surgicalproductsmag.com/product-releases/2005/03/laparoscopic-instrument-cannula-cleaning-device?qt-recent_blogs_articles=0&qt-most_popular=0