

## Patient Response Monitor

Published on Surgical Products (<http://www.surgicalproductsmag.com>)

---

## Patient Response Monitor



The B-SMART<sup>®</sup> sensory motor activated response timer is designed for use with patients undergoing conscious sedation or regional anesthesia. It evaluates neurological status by monitoring a patient's ability to detect and respond to a timed series of vibratory and auditory stimuli. Attending personnel are instantly alerted to potential problems. The B-SMART can be set for 30, 60, or 120 second repetitive time signals.

**Source URL (retrieved on 01/27/2015 - 6:14pm):**

[http://www.surgicalproductsmag.com/product-releases/2005/06/patient-response-monitor?qt-recent\\_content=0&qt-recent\\_blogs\\_articles=0](http://www.surgicalproductsmag.com/product-releases/2005/06/patient-response-monitor?qt-recent_content=0&qt-recent_blogs_articles=0)