

Decapping Pliers



The Kebby pliers decapper is designed to remove crimped aluminum seals or caps from vials and bottles. The tool is operated manually. The cap is inserted between the jaws, and removed by squeezing and tipping the pliers. This provides an easy means of removing the cap or seal safely. The stainless steel pliers feature a cushioned grip and are available in five standard sizes.

Source URL (retrieved on 01/28/2015 - 11:33am):

http://www.surgicalproductsmag.com/product-releases/2007/07/decapping-pliers?qt-recent_blogs_articles=0&qt-recent_videos=0