

Medical Suction/Vacuum Regulators



Push-To-Set[®] (PTS) suction/vacuum regulators provide a feature that automatically occludes the vacuum flow while the device is being adjusted so that the gauge accurately displays the peak vacuum. When the adjustment knob is released, the system returns to flow. The end result is the new family of PTS suction regulators will provide the lowest maximum pressure with the highest flow while potentially preventing tissue damage. In addition, the regulators incorporated a Quick-To-Max feature which requires only two turns to achieve full wall vacuum and a Stat-Set-Up feature which eliminates waiting by including an intermittent mode that starts in the "ON" cycle.

Source URL (retrieved on 02/01/2015 - 3:40pm):

http://www.surgicalproductsmag.com/product-releases/2007/09/medical-suction/vacuum-regulators?qt-recent_videos=0&qt-most_popular=0