

## OR & Procedure Basins



The new Identify-A-Bowl<sup>®</sup> OR and Procedure Basin was designed to help quickly identify solutions in the sterile fields of ORs, avoiding tragic or catastrophic mix-ups. The Identify-A-Bowl basin has a wide horizontal flange, providing a large labeling area easily seen by practitioners. The bowl has a 150-ml capacity, yet occupies minimal space on procedure trays or in sterile fields. It is sturdy so it won't tip or spill and it can be sterilized by EtO, Autoclave or Radiation.

**Source URL (retrieved on 03/06/2015 - 3:42pm):**

[http://www.surgicalproductsmag.com/product-releases/2007/09/or-procedure-basins?qt-recent\\_blogs\\_articles=0](http://www.surgicalproductsmag.com/product-releases/2007/09/or-procedure-basins?qt-recent_blogs_articles=0)