

Liquid Medical Waste Disposal System



The SAF-T Pump[®] system quickly and safely disposes of liquid medical waste directly into the sanitary sewer with no pouring required. It connects directly to existing plumbing and allows healthcare workers to dispose of liquid medical waste without exposure. Compared to disposing of full canisters, which results in more costly red bag waste, the system is also economical. The SAF-T Pump accommodates all Medi-Vac[®] suction canisters and provides a closed system that is safe, simple and economical.

Source URL (retrieved on 01/26/2015 - 2:01am):

http://www.surgicalproductsmag.com/product-releases/2007/10/liquid-medical-waste-disposal-system?qt-recent_blogs_articles=0&qt-recent_videos=0