

## Multi-Enzyme Medical Device Detergent



The KIMBERLY-CLARK\* Multi-Enzyme Detergent is available in a one-gallon pump bottle and five-gallon refills. The detergent helps assure that devices and instruments are clean before being disinfected or sterilized by penetrating and breaking down the most common contaminants found on medical devices, including proteins, fats and carbohydrate-based soils. The synergistic formulation of enzymes plus detergent works quickly within a range of water temperatures and its low-foaming property makes it appropriate for manual cleaning and for use in small washers. In addition, to help prevent waste that can result from over-measuring, the pump for the multi-enzyme detergent accurately dispenses just one ounce per pump.

**Source URL (retrieved on 02/01/2015 - 9:27pm):**

[http://www.surgicalproductsmag.com/product-releases/2008/04/multi-enzyme-medical-device-detergent?qt-recent\\_blogs\\_articles=0](http://www.surgicalproductsmag.com/product-releases/2008/04/multi-enzyme-medical-device-detergent?qt-recent_blogs_articles=0)