

Disposing Infectious Waste



Quick-Drain Systems provide safe and cost-effective disposal of liquid infectious waste. These systems satisfy OSHA, EPA, CDC and NIOSH guidelines for the removal of infectious waste risks. Two separate systems utilize Bemis Hi-Flow Rigid Suction Canisters and Bemis Quick-Fit Suction Liners. However, both can drain any competitive large-volume canister (12,000 cc and above.) Additional features include:

- The ability for waste to go into the sanitary sewer at a rate of 500 cc per second. The ability to help reduce red bag waste and eliminate per-use disposal costs such as solidifiers.
- The ability to connect to existing plumbing and not require electricity.

For more information, visit www.quick-drain.com [1]

Source URL (retrieved on 01/30/2015 - 7:52pm):

http://www.surgicalproductsmag.com/product-releases/2011/03/disposing-infectious-waste?qt-recent_content=0&qt-most_popular=0

Links:

[1] <http://www.quick-drain.com>