

Foot Control For Electrosurgical Generators



The Linemaster twin monopolar foot operated switch was designed to offer surgeons a superior accessory when performing electrosurgical procedures. The use of a foot control on electrosurgical generators combats hand fatigue that may be experienced when using surgical hand pieces. Using a foot control for soft tissue ablation applications will allow for more precise operation of the hand piece. Add-on push buttons for bipolar/monopolar combination switches as well as wireless solutions are available as required. Alternative multi-pedal or button configurations are also available to accommodate new and emerging electrosurgical functions such as tool toggle, blending and proprietary applications.

For more information, visit www.linemaster.com [1]

Source URL (retrieved on 01/26/2015 - 11:20pm):

http://www.surgicalproductsmag.com/product-releases/2011/08/foot-control-electrosurgical-generators?qt-recent_blogs_articles=0

Links:

[1] <http://www.linemaster.com/>