

High Level Disinfectant

Published on Surgical Products (<http://www.surgicalproductsmag.com>)

High Level Disinfectant



STERIS Corporation offers new Resert® XL HLD High Level Disinfectant. According to the company, the disinfectant acts as a safe alternative to aldehyde-based high level disinfectants. Features include:

- Fast-acting 8-minute high level disinfection with effective kill activity at 20°C/68°F.
- The ability to be reused for up to 21 days.
- An oxidative chemistry that eliminates staff exposure hazards associated with aldehyde-based HLD solutions (glut and OPA).
- Patented AHP® technology; 2 percent hydrogen peroxide.
- Virtually odorless, requires no special venting, measuring or mixing.
- Requires no activation or heating.
- Requires only one rinse.
- Formulated for free-rinsing, leaves no toxic residues on devices.
- The availability of Verify® chemical monitoring strips to help assure the solution is at or above the recommended concentration of hydrogen peroxide.
- Formulated with biodegradable ingredients and passes aquatic toxicity standards.
- Can be disposed of without detoxification.

High Level Disinfectant

Published on Surgical Products (<http://www.surgicalproductsmag.com>)

For more information, visit www.steris.com [1].

Source URL (retrieved on *01/30/2015 - 11:45pm*):

http://www.surgicalproductsmag.com/product-releases/2011/11/high-level-disinfectant?qt-recent_content=0&qt-recent_blogs_articles=0

Links:

[1] <http://www.steris.com/>