

Fluid Warming And Monitoring Device



Medical Solutions has introduced the 16N1, a fluid warming and monitoring device. The 16N1 provides printed documentation for the length of time that bags are warmed, and provides uniform heating of fluids. An audible and visual alarm will activate when the fluid temperature in the bag is 2°F over the set point and will continue to alarm until the bag is removed from the 16 N1 or the fluid temperature of the bag drops to within 2°F of the Set Point.

For more information, visit www.medicalsolutionsinc.com [1]

Source URL (retrieved on 01/28/2015 - 11:24pm):

http://www.surgicalproductsmag.com/product-releases/2012/06/fluid-warming-and-monitoring-device?qt-most_popular=0&qt-recent_videos=0

Links:

[1] <http://www.medicalsolutionsinc.com/>

Fluid Warming And Monitoring Device

Published on Surgical Products (<http://www.surgicalproductsmag.com>)
