

## ProneCam Monitoring System



MHP Medical's ProneCam™ Monitoring System monitors the patient's eyes and face while in the prone position. The camera cartridge slides into the front of the disposable face rest illuminating the eyes and face. Fully enclosed, the camera cartridge can be wiped clean after each procedure. The universal C-clamp allows the 4.3 inch LCD monitor to be attached to any side rail, table, IV pole or anywhere needed for patient monitoring.

For more information, visit [www.mhpmedical.com](http://www.mhpmedical.com) [1].

### Source URL (retrieved on 01/29/2015 - 12:45pm):

[http://www.surgicalproductsmag.com/product-releases/2012/09/pronecam-monitoring-system?qt-recent\\_blogs\\_articles=0&qt-recent\\_content=0](http://www.surgicalproductsmag.com/product-releases/2012/09/pronecam-monitoring-system?qt-recent_blogs_articles=0&qt-recent_content=0)

### Links:

[1] <http://www.mhpmedical.com/>