

## Decappers

Kebby Industries



Kebby Vial Decappers are used to remove crimped aluminum seals from vials in a safe manner with little effort. Seals are removed by simply gripping the seal using a tipping or rolling motion. They are available in 5 standard sizes with or without color coded latex free cushioned grips. The product is constructed of stainless steel construction with a lifetime guarantee.

For more information, visit

[https://www.vialcrimper.com/Plier\\_Decapper\\_decrimping.html](https://www.vialcrimper.com/Plier_Decapper_decrimping.html) [1] and [www.kebbyindustries.com](http://www.kebbyindustries.com) [2]

**Source URL (retrieved on 01/29/2015 - 4:50pm):**

[http://www.surgicalproductsmag.com/product-releases/2013/04/decappers?qt-recent\\_videos=0&qt-recent\\_content=0](http://www.surgicalproductsmag.com/product-releases/2013/04/decappers?qt-recent_videos=0&qt-recent_content=0)

### Links:

[1] [https://www.vialcrimper.com/Plier\\_Decapper\\_decrimping.html](https://www.vialcrimper.com/Plier_Decapper_decrimping.html)

[2] <http://www.kebbyindustries.com/>