

Single-Use Active Self-Warming Blanket

Mölnlycke Health Care



The BARRIER EasyWarm from Mölnlycke Health Care active self-warming blanket has 12 air-activated warming pads consisting of 100 percent natural material such as active coal, clay, salt, water, and iron powder. The blanket maintains an average temperature of 44 degrees Celsius for up to 10 hours, raising skin temperature to a maximum of 42 degrees Celsius. Quick and easy to set up, the blanket can be used before, during and after surgery as it does not need any additional equipment such as fans or hoses of electricity.

For more information, visit www.molnlycke.com [1]

Source URL (retrieved on 01/28/2015 - 9:14am):

<http://www.surgicalproductsmag.com/product-releases/2013/12/single-use-active-self-warming-blanket>

Links:

[1] <http://www.molnlycke.com/>